

PowerPoint / Keynote and OsiriX Integration
Dr Sameer Shamshuddin, Royal Lancaster Infirmary,
United Kingdom

How to Download, Install and Prepare your Talk:

1. Visit Turyon website and download the software:
<http://www.turyon.com/downloads/SimpleHTTPServer.app.zip>
2. Copy downloaded zip file to Applications and unzip the contents
3. Launch from applications-SimpleHTTPServer App
4. Enter the Osirix Port Number on the right hand side---(8080)enter
(Tip- You can find the osirix port number from Osirix-Pref-Listner)
5. Quit application
6. Launch keynote and safari
7. In safari / Chrome type/ paste the link:
<http://localhost:12345/keynote?type=patientID&name=SAMSHA-0280&close=1>
8. Change SAMSHA-0280 to your case code (PatientID in your database)
9. Save your link URL as a favorite (will help copy paste to Keynote in future)
10. In Keynote, you can highlight any text and add a link through inspector (type/copy the link above, but change the PatientID to the case you intend to show and make sure that the actually case exist in the opened database)
11. In SlideShow mode, just click the link to launch your desired case from osirix

Conference Settings

1. Launch OsiriX, Keynote, Chrome/Safari and then start the ServerApp
2. In Keynote or PowerPoint, go to slideshow mode
3. Click the link in the presentation and it will launch the case within OsiriX

Important - Make sure to keep the cases in the default database or keep the database containing the case open, if you have multiple databases

4. Switch back to presentation with **CMD+Tab**

Note from the developers:

Please consider making a donation by Paypal (turyon@turyon.com)

Thank you for your interest, I hope you'll count on us for future projects and developments.

Juan Jose Cermeño Portet

Turyon Projects S.L, CTO
Medical Imaging Consultant, Software Developer
email: jjcermeno@turyon.com
www.turyon.com

Happy days ahead, stun your audience 😊

Sameer Shamshuddin 😊